

JEOL 100CX II Startup and Alignment

Always

- Change samples or film with filament current **OFF**.
- Handle any parts that go into the column (film, film cassettes, film boxes, specimens, specimen holder) with gloved hands or tweezers.
- Call for help if you've lost the beam and you're not sure how to find it.
- Be gentle. **No** force is required to operate this TEM.

Inserting the Specimen

- Insert 3mm specimen into **Specimen Holder Rod**. **Gently** tap the holder near the specimen to demonstrate that it is held securely. Inspect holder o-ring for dirt
– clean and regrease if nec.
- Set specimen position to 1 or 2 as appropriate.
- Insert specimen holder into goniometer, matching pin on holder to notch at 9 o'clock. Hold in with slight pressure until red vacuum indicator lights up. When light goes out, turn holder clockwise and guide it into position. **Do not release the holder until it is fully in position.**

Startup Procedure

- If the **Ready**, **DP** or **High Vac** lights are off, call for help.
- If necessary, fill the **Cold Trap** with liquid nitrogen.
- Make sure that the **Objective Aperture** and **SAD Aperture** are out.
- Select the **Mag Function**.
- Select a **Magnification** of ~2 kx.
- Bring up the High Tension (HT):
 - If the HT is off (**HT** button is unlit):
 - Select an **Accelerating Voltage** of 20 kV.
 - Push in the **HT** button.
 - When the **Beam Current** indicator stabilizes, (this may take a long time) select 40 kV.
 - Continue ramping up to the desired Accelerating Voltage. You may need to jump from 40kV directly to 100kV
 - If the **HT** is on, but at an accelerating voltage higher than desired:
 - Push in the desired **HT**. You need not ramp downwards.
 - If the **HT** is on, but at an accelerating voltage lower than desired:
 - Ramp the voltage up one step (+20kV) at a time, waiting at least 1 minute for the **Beam Current** meter to stabilize before each increase in HT. You may need to skip the 60 and 80 positions if the emission will not stabilize.
- Note the Beam Current reading (dark current). Turn the **Filament Emission** knob clockwise until the pointer meets the metal stop. **Do not adjust the position of the stopper.** If the Beam Current increases more than 20 ma above dark current, call for help

- **If you don't see the beam:**
 - Fully spread the beam by turning the **Condenser** knob clockwise.
 - Double check that the **Objective** and **SAD Apertures** are out
 - Switch to the **Low Mag** Function.
 - Using the map on the upper left panel, systematically search for the specimen perforation.
 - Remove the **Specimen Holder Rod**, or select a holder position that has no specimen
 - Call for help

Routine Alignment

Gun Alignment With **Magnification** at ~5 kx, **Spot Size** at 1, **Function** at **Mag**, bring **Condenser** to cross-over. Desaturate the **Filament** and center the filament image with the **Gun Align Tilt** knobs (lower right panel). Resaturate **Filament**.

Condenser Alignment With **Magnification** at ~5 kx, **Spot Size** at 3, bring **Condenser** to cross-over. Center beam with **Alignment Trans** knobs (upper right panel). Switch **Spot Size** to 1, bring **Condenser** to cross-over, center beam with **Gun Align Trans knobs** (lower right panel). Repeat until spots 1 and 3 coincide.

Condenser Aperture Centering With **Magnification** at ~5 kx, insert desired **Condenser Aperture**. Bring **Condenser** to cross-over and center beam with **Alignment Trans** knobs. With **Condenser**, spread beam clockwise and center beam with **Condenser Aperture x** and **y** controls. Repeat until beam spreads concentrically.

Condenser Stigmation With **Magnification** at ~10 kx, bring **Condenser** to cross-over. Desaturate **Filament**. Adjust **Condenser Stigmator** knobs (lower right panel) until you obtain a sharp filament image. Resaturate **Filament**.

Z Axis (height) Correction With **Magnification** at ~10 kx, center a recognizable specimen feature. Rock the **Goniometer** back and forth. Adjust **Z Adjust** knob until the image motion is minimized. Refocus the image.

Deflector Coil Balancing With **Magnification** at ~10 kx, bring **Condenser** to cross-over. On lower left panel, switch **Wobbler Up** to **X**, then converge the 2 spots with **X Comp** and **X Corr** knobs. Repeat for **Y**.

Image Wobbler Balancing With the image focused, bring **Condenser** to cross-over and turn on the **Image Wobbler** switch (on upper left panel). Converge the 2 spots with the **Image Wobbler A** and **B** knobs (on lower right panel).

Objective Lens Alignment (Voltage Centering) With **Magnification** > 30,000x, beam centered, center a recognizable specimen feature. Switch on **HV Wobbler** (on upper left panel) and engage the focusing binoculars. Minimize image movement with the **Align Tilt** knobs (on upper right and left). Recenter illumination as necessary with **Alignment Trans** knobs upper right and left).

Objective Aperture Centering Insert an **SAD Aperture** and switch **Function** to **SA Diff**. Insert the desired **Objective Aperture** and focus the aperture image. Center the aperture around the forward scattered spot. Return **Function** to **Mag**.

Objective Lens Stigmation Select **Obj Stigmator** set #1 or #2 with the switch on the lower left panel. For rough stigmation, underfocus or overfocus the image, then balance the Fresnel fringes around a specimen hole with **Coarse** and **Fine X** and **Y** controls (on lower left). Above 200 kx, find **Focus** of minimum contrast. Further minimize contrast with the **Fine X** and **Y** controls.